Specifications

System		Value/Data
	RF Power	Up to 50W
	RF Frequency	1MHz

TH TOWER	Op to sow
RF Frequency	1MHz
Input voltage and current	100-240 Volt, 50-60 Hz, max 2.2A
User interface	LCD color touch screen
Dimensions	L45 x W45 x H110 cm; L17.7 x W17.7 x H43.3 in.
Weight	~30 Kg; ~66 Lb

Applicator

	Hybrid Energy [™] applicator	TriLipo® RF applicators
Number of needles/tip	Gen 36: 6x6 Gen 12: 2x6	n/a
Disposable tip pulses	Up to 800 pulses per tip	n/a
Number of electrodes	n/a	Large – 5, Medium – 3, Small - 3

As Featured in the Media



About Lumenis

Lumenis is a global leader in the field of minimally-invasive clinical solutions for the Surgical, Ophthalmology and Aesthetic markets, and is a world-renowned expert in developing and commercializing innovative energy-based technologies, including Laser, Intense Pulsed Light (IPL) and Radio-Frequency (RF). For nearly 50 years, Lumenis' ground-breaking products have redefined medical treatments and have set numerous technological and clinical gold-standards. Lumenis has successfully created solutions for previously untreatable conditions, as well as designed advanced technologies that have revolutionized existing treatment methods in each and every one of the verticals we operate in. Our drive for innovation stems from an uncompromising commitment to improving the health and well-being of our patients; addressing new and growing needs of aging populations; and in offering medical professionals cutting-edge solutions that fit seamlessly into the health-economics environment of the 21st century. The world over, we bring Energy to Healthcare.

About Pollogen

Pollogen Ltd. is a global leader in the medical aesthetics market providing innovative, safe, and effective solutions for women and men that enhance and expand the practice of medical aesthetics professionals. The company offers products to address a range of skin treatments under the industry's six premier brands: TriPollar®, TriLipo®, TriFractional™, OxyGeneo®, Hybrid Energy™, and VoluDerm™


Pollogen offers a full line of clinically-proven, non-invasive anti-aging facial and body contouring treatment platforms for circumference reduction, cellulite reduction and skin tightening. The company's aesthetic devices: Pollogen LEGEND, geneO+, Maximus, Apollo, and Surgen offer customers clinically proven, safe and effective anti-aging lunchtime beauty treatments that deliver immediate and long-lasting results. With a distribution network in over 60 countries, Pollogen is providing innovative solutions to medical professionals globally.

© 2016 Pollogen Ltd. All rights reserved. Trademarks and logos of Apollo, RegenXL, TopJen, Maximus, geneO+, Pollogen LEGEND, NeoRevive, NeoBright, Surgen, TriPollar, TriLipo, TriFractional, VoluDerm, Hybrid Energy and Pollogen are all trademarks of Pollogen Ltd. TriPollar, TriLipo, TriFractional, VoluDerm and






Hybrid Energy™ with volumizing effect



RF energy is applied, providing a safe and effective penetration for controlled heating targeting deep dermis.



2nd electrical energy enhances the total clinical effect - directly applied to the specific target area creates microwounds triggering production of new collagen, elastin & hyaluronic acid.

- Dermal volumizing
- · Skin tightening and fine line reduction
- · Wrinkle reduction
- · Skin resurfacing
- · Facial skin contouring

TriLipo® RF energy + Dynamic Muscle Activation (DMA)



RF Energy heats the subcutaneous fat tissue layer stimulating the release of liquid fat from the cells accelerating natural fat metabolism.



Simultaneous internal and external pressure on the fat layer drains the fat via the DMA lymphatic drainage mechanism optimizing tissue oxygenation and detoxification.

- · Wrinkle reduction
- · Cellulite and circumference reduction
- · Body contouring
- · Skin firming, lifting & rejuvenation

A WINNING Combination



MEET Your Patients' Needs

Hybrid Energy[™] with volumizing effect

- 1. Fixes those hard to treat areas: eyes*, upper lip and neck
- 2. Wrinkle reduction & dermal volumizing after only one treatment
- 3. A single treatment takes 30 minutes. Ideal for the busy woman

TriLipo® with DMA

- 1. Non-invasive body sculpting & skin tightening
- 2. Clinically proven & effective for all skin types
- 3. Walk-in /Walk-out procedure

The Pollogen LEGEND provides non-surgical solutions that simply work. With each additional treatment, results are enhanced, there is no down time and the treatment is painless.

Face & Body



VISIBLE Immediate Long-term RESULTS

Hybrid Energy with volumizing effect









Volumetric increase & wrinkle reduction Courtesy of Steven D. Shapiro M.D., Palm Beach Gardens, FL.

Dermal volumizing and wrinkle reduction Courtesy of Dr. Alex Levenberg, M.D., Tel Aviv, Israel

TriLipo with DMA











After 5 treatments

Thighs cellulite reduction Courtesy of Dvora Ancona, M.D., Milan, Italy

Abdomen skin tightening Courtesy of Dr. Alex Levenberg, M.D., Tel Aviv, Israel



"With the Pollogen LEGEND+ system, I can totally achieve the results my patients and I desire. The facial dermal volumizing and wrinkle reduction results are exceptional. The skin becomes more supple and elastic and my patients look and feel years younger. Most importantly, the results are natural-looking which is what my patients want."

Michael H. Gold, MD, Nashville, TN



"I find Pollogen's TriLipo non-invasive body shaping systems effective for circumference reduction, removing localized fat deposits and reducing cellulite. My patients are very pleased with the immediate and long-lasting results and are grateful that the treatments are so quick and pleasant."

Michele Pelletier, M.D. Toulon France

Clinical Validation

Wanitphakdeedecha R. Treatment of abdominal cellulite and circumference reduction with radiofrequency and dynamic muscle activation. Journal of Cosmetic and Laser Therapy. 2015; 17:5

"In conclusion, circumference and cellulite can be significantly and safely reduced with the use of the RF and dynamic muscle activation technology."

Shapiro, S. VoluDerm Micro-Needle Technology for Treating Skin Laxity and Wrinkles-Initial Clinical Experience. Jacobs Journal of Experimental Dermatology, Jan 2015

"VO clinical results demonstrate skin texture improvement, volumizing effects and improved wrinkle appearance."

Gershonowitz A, Gat A. VoluDerm Microneedle Technology for Skin Treatments - In Vivo Histological Evidence. Journal of Cosmetic and Laser Therapy, Posted online on September 2014

"It was concluded that exposure to the HE technology was shown to affect significant collagen remodeling and a significant increase in mid- to deep dermis GAGs..."

Boisnic S, Branchet MC. Ex-vivo Study of Hybrid Energy Technology using a Human Skin Model. European Journal of Dermatology. 2014 Jan-Feb; 24(1):46-52

"Results demonstrated a significant increase of the epidermal mitotic index, significant dermal collagen remodeling and significant increases in both epidermal and dermal GAGs."

Potekaev, N, Zhukova, O. Evaluation of Safety and Efficacy of the Maximus™ System for Facial Wrinkles. Journal of Cosmetics, Dermatological Sciences and Applications, 2013, 3, 151-156

"Conclusions: The Maximus system powered by the TriLipo technology is a non-invasive, effective, safe, and virtually painless treatment for reduction of wrinkles and facial skin rejuvenation."

Manuskiatti W. TriPollar Radiofrequency. Cellulite, Pathophysiology and Treatment second edition. Goldman MP, Hexsel D., 2010; 22:158-167

